



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EnzyChrom(TM) Fumarate Assay Kit BAS-EFUM-100

Article Name	EnzyChrom(TM) Fumarate Assay Kit
Biozol Catalog Number	BAS-EFUM-100
Supplier Catalog Number	EFUM-100
Alternative Catalog Number	BAS-EFUM-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of fumarate and drug effects on its metabolism. Key Features: Fast and sensitive. Use of 20 μ L sample. Linear detection range 0.005 to 2 mM fumarate in 96-well plate assay. Convenient. The procedure involves adding a si