



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) Galactose Assay Kit BAS-EGAL-100

Article Name	EnzyChrom(TM) Galactose Assay Kit
Biozol Catalog Number	BAS-EGAL-100
Supplier Catalog Number	EGAL-100
Alternative Catalog Number	BAS-EGAL-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of galactose and evaluation of drug effects on its metabolism. Key Features: Use as little as 20 $\mu L$ samples. Linear detection range in 96-well plate: 10 to 1000 $\mu M$ galactose for colorimetric assays and 10 to 100 $\mu M$ for