



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) Acetate Assay Kit BAS-EOAC-100

Article Name	EnzyChrom(TM) Acetate Assay Kit
Biozol Catalog Number	BAS-EOAC-100
Supplier Catalog Number	EOAC-100
Alternative Catalog Number	BAS-EOAC-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of acetic acid or acetate and evaluation of drug effects on acetate metabolism. Key Features: Sensitive. Use as little as 10 $\mu$ L samples. Detection range: 0.20 to 20 mM acetate for colorimetric assays and 0.13 to 2 mM fo