



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EnzyChrom(TM) Sialic Acid Assay Kit BAS-ESLA-100

Article Name	EnzyChrom(TM) Sialic Acid Assay Kit
Biozol Catalog Number	BAS-ESLA-100
Supplier Catalog Number	ESLA-100
Alternative Catalog Number	BAS-ESLA-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of free and total sialic acid (NANA) and evaluation of drug effects on its metabolism. Key Features: Sensitive and accurate. Use as little as 10 μ L samples. Linear detection range in 96-well plate: 0.02 to 1 mM sialic a