



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## EnzyChrom(TM) Xanthine Assay Kit BAS-EXAN-100

Article Name	EnzyChrom(TM) Xanthine Assay Kit
Biozol Catalog Number	BAS-EXAN-100
Supplier Catalog Number	EXAN-100
Alternative Catalog Number	BAS-EXAN-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For quantitative determination of xanthine and drug effects on its metabolism. Key Features: Sensitive and accurate. Use as little as 10 $\mu$ L samples. Linear detection range in 96-well plate for 30 minute incubation: 0.01 to 2.5 mM xanthine for colorim