



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## QuantiChrom(TM) Glucose-6-Phosphate Dehydrogenase Inhibitor Screening Kit BAS-IGPDH-100

Article Name	QuantiChrom(TM) Glucose-6-Phosphate Dehydrogenase Inhibitor Screening Kit
Biozol Catalog Number	BAS-IGPDH-100
Supplier Catalog Number	IGPDH-100
Alternative Catalog Number	BAS-IGPDH-100
Manufacturer	BioAssay Systems
Category	Kits/Assays
Product Description	For evaluation and high-throughput screen (HTS) of glucose-6-phosphate dehydrogenase modulators. Key Features: Safe. Non-radioactive assay. High-throughput. Homogenous "mix-incubate-measure type assay. Can be readily automated on HTS liquid handling