



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

PhosphoBlocker(TM) Blocking Reagent CBL-AKR-103

Article Name	PhosphoBlocker(TM) Blocking Reagent
Biozol Catalog Number	CBL-AKR-103
Supplier Catalog Number	AKR-103
Alternative Catalog Number	CBL-AKR-103
Manufacturer	Cell Biolabs
Category	Sonstiges
Product Description	Most commercially available Western blot blockers, such as dry milk or serum, are sufficient to block unreactive sites on the membrane. However, they are not designed to preserve phosphoprotein antigens during blotting. Our PhosphoBLOCKER(TM) Blockin