



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **Total Phosphatidic Acid Assay Kit (Fluorometric) CBL-MET-5019**

Article Name	Total Phosphatidic Acid Assay Kit (Fluorometric)
Biozol Catalog Number	CBL-MET-5019
Supplier Catalog Number	MET-5019
Alternative Catalog Number	CBL-MET-5019
Manufacturer	Cell Biolabs
Category	Kits/Assays
Product Description	Our Total Phosphatidic Acid Assay Kit measures total phosphatidic acid content (PA and LPA) in cell lysate samples by a coupled enzymatic reaction system. First, lipase is used to hydrolyze phosphatidic acid to glycerol-3-phosphate. Next, the glycero