



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

DAG (Diacylglycerol) Assay Kit CBL-MET-5028

Article Name	DAG (Diacylglycerol) Assay Kit
Biozol Catalog Number	CBL-MET-5028
Supplier Catalog Number	MET-5028
Alternative Catalog Number	CBL-MET-5028
Manufacturer	Cell Biolabs
Category	Kits/Assays
Product Description	Our DAG (Diacylglycerol) Assay Kit measures diacylglycerol in cell lysate samples by a coupled enzymatic reaction system. First, kinase is used to phosphorylate DAG in samples, yielding phosphatidic acid. Next, lipase is used to hydrolyze phosphatidi