



Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

OxiSelect(TM) ORAC Activity Assay Kit CBL-STA-345

Article Name	OxiSelect(TM) ORAC Activity Assay Kit
Biozol Catalog Number	CBL-STA-345
Supplier Catalog Number	STA-345
Alternative Catalog Number	CBL-STA-345
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	Fluorometric
Product Description	The ORAC assay is a powerful tool to measure the antioxidant capacity of biomolecules. Our OxiSelect(TM) ORAC Activity Assay measures this capacity in less than 2 hours from a wide variety of sample types