

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## OxiSelect(TM) Human Oxidized LDL ELISA Kit (CML-LDL Quantitation) CBL-STA-388

Article Name	OxiSelect(TM) Human Oxidized LDL ELISA Kit (CML-LDL Quantitation)
Biozol Catalog Number	CBL-STA-388
Supplier Catalog Number	STA-388
Alternative Catalog Number	CBL-STA-388
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	Colorimetric
Product Description	LDL, so-called bad cholesterol, is even more dangerous when it becomes oxidized. Oxidized LDL (OxLDL) is more reactive with surrounding tissues and can collect within the inner-lining of arteries. Macrophages, cholesterol, and other lipids can accumu