

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

OxiSelect(TM) Human Oxidized LDL ELISA Kit (HNE-LDL Quantitation) CBL-STA-389

Article Name	OxiSelect(TM) Human Oxidized LDL ELISA Kit (HNE-LDL Quantitation)
Biozol Catalog Number	CBL-STA-389
Supplier Catalog Number	STA-389
Alternative Catalog Number	CBL-STA-389
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	Colorimetric
Product Description	LDL, so-called bad cholesterol, is even more dangerous when it becomes oxidized. Oxidized LDL (OxLDL) is more reactive with surrounding tissues and can collect within the inner-lining of arteries. Macrophages, cholesterol, and other lipids can accumu