

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## QuickTiter(TM) Lentivirus Titer Kit (Lentivirus-Associated p24 ELISA) CBL-VPK-107

Article Name	QuickTiter(TM) Lentivirus Titer Kit (Lentivirus-Associated p24 ELISA)
Biozol Catalog Number	CBL-VPK-107
Supplier Catalog Number	VPK-107
Alternative Catalog Number	CBL-VPK-107
Manufacturer	Cell Biolabs
Category	Kits/Assays
Application	Colorimetric
Product Description	Measuring the HIV-1 p24 antigen is a long-established method for lentivirus quantitation. However, the traditional p24 ELISA detects both virus-associated p24 and free p24 generated by 293 cells during transient transfection. Free p24 can account for