

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

QuickTiter(TM) Lentivirus Titer Kit (Lentivirus-Associated HIV p24), Trial Size CBL-VPK-107-T

| Article Name | QuickTiter(TM) Lentivirus Titer Kit (Lentivirus-Associated HIV p24), Trial Size |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Biozol Catalog Number | CBL-VPK-107-T |
| Supplier Catalog Number | VPK-107-T |
| Alternative Catalog Number | CBL-VPK-107-T |
| Manufacturer | Cell Biolabs |
| Category | Kits/Assays |
| Application | Colorimetric |
| Product Description | Measuring the HIV-1 p24 antigen is a long-established method for lentivirus quantitation. However, the traditional p24 ELISA detects both virus-associated p24 and free p24 generated by 293 cells during transient transfection. Free p24 can account for |