

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Horse normal serum () AGR-AS10-1562

| | |
|----------------------------|---|
| Article Name | Horse normal serum () |
| Biozol Catalog Number | AGR-AS10-1562 |
| Supplier Catalog Number | AS10-1562 |
| Alternative Catalog Number | AGR-AS10-1562 |
| Manufacturer | Agrisera |
| Host | Equine |
| Category | Zellen/Zellkultur |
| Product Description | Normal horse serum lipid extracted and dialyzed against 10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2. This product is used as a blocking reagent or control for most immunoassay applications. Blocking is a critical step in most immunoassays ... |
| Purity | Serum |
| Form | Lyophilized |
| Antibody Type | Secondary Antibody |
| Application Notes | Usually 5 % (v/v) serum blocking buffer it is recommended and for that you need to mix 2 ml of blocking serum with 38 ml of diluent buffer (PBS with Tween 20 detergent) to obtain a total volume of 40 ml, Use immediately or store at 2-8 C or colder, |