

Diagnostica Vertrieb GmbH, Leipziger Straße 4 85386 Eching, Germany Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99 E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Horse normal serum () AGR-AS10-1562

Article Name	Horse normal serum ()
Biozol Catalog Number	AGR-AS10-1562
Supplier Catalog Number	AS10-1562
Alternative Catalog Number	AGR-AS10-1562
Manufacturer	Agrisera
Host	Equine
Category	Zellen/Zellkultur
Product Description	Normal horse serum lipid extracted and dialyzed against 10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2. This product is used as a blocking reagent or control for most immunoassay applications. Blocking is a critical step in most immunoassays
Purity	Serum
Form	Lyophilized
Antibody Type	Secondary Antibody
Application Notes	Usually 5 % (v/v) serum blocking buffer it is recommended and for that you need to mix 2 ml of blocking serum with 38 ml of diluent buffer (PBS with Tween 20 detergent) to obtain a total volume of 40 ml, Use immediately or store at 2-8 C or colder,