

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit F(ab)2 anti-Guinea Pig IgG (F(ab)2)-unconj., MinX none, Polyclonal DNA-SEC-182672

|                            |   |
|----------------------------|---|
| Article Name               | Rabbit F(ab)2 anti-Guinea Pig IgG (F(ab)2)-unconj., MinX none, Polyclonal   |
| Biozol Catalog Number      | DNA-SEC-182672  |
| Supplier Catalog Number    | SEC-182672  |
| Alternative Catalog Number | DNA-SEC-182672  |
| Manufacturer               | dianova   |
| Host                       | Rabbit  |
| Category                   | Antikörper  |
| Application                | WB, IHC, ELISA  |
| Species Reactivity         | Guinea pig  |
| Immunogen                  | Guinea Pig IgG F(ab)2 fragment  |
| Conjugation                | Unconjugated  |
| Product Description        | F(ab)2 Anti-Guinea Pig IgG F(ab)2 Antibody generated in rabbit detects Guinea Pig F(ab)2. Representing approximately 75% of serum immunoglobulins, IgG is the most abundant antibody isotype found in the circulation. IgG molecules are synthesized and s... |
| Clonality                  | Polyclonal  |
| Isotype                    | Ig  |
| Buffer                     | 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2  |

|                    |   |
|--------------------|---|
| Purity             | This product is a F(ab)2 fragment of IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by chromatographic separations |
| Formula            | 20 mM K3PO4, 150 mM NaCl, pH 7.2, lyophilisate, Azide/BSA free  |
| Target             | Guinea Pig  |
| Antibody Type      | Polyclonal Antibody   |
| Application Dilute | WB: 1:2,000-1:10,000  |
| Application Notes  | Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring extremely low background levels, absence of Fc mediated binding, lot-to-lot consistency, high titer and specificity. The maxi      |