

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit F(ab)2 anti-Sheep IgG (H+L)-HRPO, MinX none, Polyclonal DNA-SEC-182717

Article Name	Rabbit F(ab)2 anti-Sheep IgG (H+L)-HRPO, MinX none, Polyclonal
Biozol Catalog Number	DNA-SEC-182717
Supplier Catalog Number	SEC-182717
Alternative Catalog Number	DNA-SEC-182717
Manufacturer	dianova
Host	Rabbit
Category	Antikörper
Application	WB, IHC, ELISA
Species Reactivity	Sheep
Immunogen	Sheep IgG whole molecule
Conjugation	HRP
Product Description	F(ab)2 Anti-Sheep IgG (H&L) Peroxidase Antibody generated in rabbit detects immunoglobulin g from sheep, both heavy and light chains of the antibody molecule are present. Each IgG has two antigen binding sites. Representing approximately 75% of serum
Clonality	Polyclonal
Concentration	10.0 mg/mL
Isotype	Ig
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Purity	This product is a F(ab)2 fragment of an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by extensive dialysis a
Formula	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% Gentamicin
Target	Sheep
Antibody Type	Polyclonal Antibody
Application Dilute	ELISA Dilution: 1:10,000 - 1:50,000, Immunohistochemistry Dilution: 1:500 - 1:2,500, Western Blot Dilution: 1:1,000 - 1:10,000
Application Notes	F(ab)2 Anti-Sheep IgG (H&L) Peroxidase secondary antibody reagents are ideal for ELISA, western blotting, Immunohistochemistry, Fluorescence Microscopy, Flow Cytometry as well as other antibody detection methods.