

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit IgG anti-Goat IgG (H+L)-ATTO 425, MinX Ms,Rb, ATTO-425, Polyclonal DNA-SEC-182860

Article Name	Rabbit IgG anti-Goat IgG (H+L)-ATTO 425, MinX Ms,Rb, ATTO-425, Polyclonal
Biozol Catalog Number	DNA-SEC-182860
Supplier Catalog Number	SEC-182860
Alternative Catalog Number	DNA-SEC-182860
Manufacturer	dianova
Host	Rabbit
Category	Antikörper
Application	WB, IF, FLISA
Species Reactivity	Goat
Immunogen	Goat IgG whole molecule
Conjugation	ATTO-425
Product Description	Anti-Goat IgG ATTO dye Antibody generated in rabbit detects goat IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as well
Clonality	Polyclonal
Concentration	1.0 mg/mL
Isotype	Ig

Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Purity	Goat IgG (H&L) Antibody conjugated to ATTO425 was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelec
Formula	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target	Goat
Antibody Type	Polyclonal Antibody
Application Dilute	FLISA Dilution: >1:20,000, Fluorochrome Protein Value: 2.0, IF Microscopy Dilution: >1:5,000, Western Blot Dilution: >1:10,000
Application Notes	Rabbit Anti-Goat IgG (H&L) Antibody conjugated to ATTO425 is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis