

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat IgG anti-Rabbit IgG (H+L)-FITC, MinX none, Polyclonal DNA-SEC-183357

Article Name	Goat IgG anti-Rabbit IgG (H+L)-FITC, MinX none, Polyclonal
Biozol Catalog Number	DNA-SEC-183357
Supplier Catalog Number	SEC-183357
Alternative Catalog Number	DNA-SEC-183357
Manufacturer	dianova
Host	Goat
Category	Antikörper
Application	DOT, WB
Species Reactivity	Rabbit
Immunogen	Rabbit IgG whole molecule
Conjugation	FITC
Product Description	Anti-Rabbit IgG Antibody Fluorescein generated in goat detects rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as
Clonality	Polyclonal
Concentration	2.0 mg/mL
Isotype	Ig
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Purity	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a singl
Formula	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% NaN3
Target	Rabbit
Antibody Type	Polyclonal Antibody
Application Dilute	FLISA Dilution: 1:10,000 - 1:50,000, Flow Cytometry Dilution: 1:500 - 1:2,500, Fluorochrome Protein Value: 3.97, IF Microscopy Dilution: 1:1,000 - 1:5,000
Application Notes	Anti-Rabbit IgG Antibody Fluorescein has been tested by dot blot and western blot and is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex ana