

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Goat IgG anti-Rabbit IgG (H+L)-ATTO 532, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-532, Polyclonal DNA-SEC-183383

Article Name	Goat IgG anti-Rabbit IgG (H+L)-ATTO 532, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-532, Polyclonal
Biozol Catalog Number	DNA-SEC-183383
Supplier Catalog Number	SEC-183383
Alternative Catalog Number	DNA-SEC-183383
Manufacturer	dianova
Host	Goat
Category	Antikörper
Application	DOT, WB
Species Reactivity	Rabbit
Immunogen	Rabbit IgG whole molecule
Conjugation	ATTO-532
Product Description	Anti-Rabbit IgG (H&L) ATTO 532 Antibody generated in goat detects reactivity to Rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to virus...
Clonality	Polyclonal
Concentration	1.0 mg/mL
Isotype	Ig

Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Purity	Rabbit IgG (H&L) Antibody ATTO 532 was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis
Formula	20 mM K ₃ PO ₄ , 150 mM NaCl, pH 7.2, lyophilisate, 0.01% NaN ₃
Target	Rabbit
Antibody Type	Polyclonal Antibody
Application Dilute	FLISA Dilution: >1:20,000, Fluorochrome Protein Value: 1.4, IF Microscopy Dilution: >1:5,000, Western Blot Dilution: >1:10,000
Application Notes	Anti-Rabbit IgG (H&L) conjugated to ATTO 532 has been tested by dot blot and western blot and is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is