

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Goat IgG anti-Rabbit IgG (H+L)-ATTO 550, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-550, Polyclonal DNA-SEC-183386**

Article Name	Goat IgG anti-Rabbit IgG (H+L)-ATTO 550, MinX Bo,Ck,Go,Gp,Hm,Ho,Hu,Ms,Rt,Sh, ATTO-550, Polyclonal
Biozol Catalog Number	DNA-SEC-183386
Supplier Catalog Number	SEC-183386
Alternative Catalog Number	DNA-SEC-183386
Manufacturer	dianova
Host	Goat
Category	Antikörper
Application	DOT, WB
Species Reactivity	Rabbit
Immunogen	Rabbit IgG whole molecule
Conjugation	ATTO-550
Product Description	Anti-Rabbit IgG (H&L) ATTO 550 Antibody generated in goat detects reactivity to Rabbit IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to virus...
Clonality	Polyclonal
Concentration	1.0 mg/mL
Isotype	Ig

Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Purity	Rabbit IgG (H&L) Antibody ATTO 550 was prepared from monospecific antiserum by immunoaffinity chromatography using Rabbit IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis
Formula	20 mM K <sub>3</sub> PO <sub>4</sub> , 150 mM NaCl, pH 7.2, lyophilisate, 0.01% NaN <sub>3</sub>
Target	Rabbit
Antibody Type	Polyclonal Antibody
Application Dilute	FLISA Dilution: >1:20,000, Fluorochrome Protein Value: 1.9, IF Microscopy Dilution: >1:5,000, Western Blot Dilution: >1:10,000
Application Notes	Anti-Rabbit IgG (H&L) conjugated to ATTO 550 has been tested by dot blot and western blot and is designed for STED microscopy, FRET, immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is