

Diagnostica Vertrieb GmbH, Leipziger Straße 4

85386 Eching, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Phosphatase Inhibitor Cocktail III (2 Tubes, 100x) TGM-C0004

Article Name	Phosphatase Inhibitor Cocktail III (2 Tubes, 100x)
Biozol Catalog Number	TGM-C0004
Supplier Catalog Number	C0004
Alternative Catalog Number	TGM-C0004-1MLI+1MLII,TGM-C0004-(1MLI+1MLII)X10,TGM-C0004-(1MLI+1MLII)X50,TGM-C0004-(1MLI+1MLII)X100
Manufacturer	TargetMol
Category	Biochemikalien
Product Description	It is easy for the proteins to be degraded or dephosphorylated during in vitro extraction procedure, which could cause an inaccurate result of protein expression detection. Therefore, adding protease inhibitors or phosphatase inhibitors into the extr